



Disposable Respiratory Mask - FFP2 NR MUSK001

These particulate respirators protect against solid and non-volatile liquid particles.

WARNINGS AND LIMITATIONS

Failure to follow all instructions and limitations on the use of these respirators and/or failure to wear these respirators during all times of exposure can reduce effectiveness and may result in sickness or death. For additional information please contact us

LIMITATIONS

1. These respirators meet the requirements of EN 149:2001+A1:2009, filtering face piece respirators for use against particles. They should be used to protect the wearer from solid and non-volatile liquid particles only.
2. Proper selection, training, use and appropriate maintenance are essential for the respirator in order to protect the wearer from certain airborne contaminants. Failure to follow all instructions on the use of these respiratory protection and/or failure to properly wear the complete respirator during all periods of exposure may adversely affect the wearer's health, easily leading to severe or life threatening illness or permanent disability..
3. Always check and make sure that the respirator is within the stated shelf life. (Use by Date)
4. Before use, the wearer must be trained in the use of complete respirator in accordance with health and safety standards/guidance.
5. These respirators do not contain components made from natural rubber latex.
6. These respirators do not protect against gases/vapors.
7. These respirators are not suitable for children.
8. Do not use in atmosphere containing less than 19.5% oxygen.
9. Do not use for respiratory protection against atmospheric contaminants/ concentrations which are unknown or immediately dangerous to life and health (IDLH)
10. Do not use with beards or other facial hair, which might inhibit contact between the face and the respirator ; thus preventing a good seal.
11. Do not use in the presence of high intensive heat source or flammable gas.
12. Leave the contaminated area immediately if :
 - a) Breathing becomes difficult.
 - b) Dizziness or other distress occurs.
13. Always be sure you are using the complete respirator.
14. The respirator is worn during period of exposure of intend only.
15. Replace respirator when necessary.
16. Discard and replace the respirator if it becomes damaged, breathing resistance becomes excessive or at the end of the shift.



MUSK 001

[STANDARDS]

EN149 : 2001 + A1 : 2009

REGULATION(EU)2016/425

ISO 13485 : 2016

USUAGE METHOD



1. Clean your hands with soap and water before touching the mask.



2. With noseclip facing up, unfold the mask with both hands and pull the elastic over the ears.



3. Mold the nose strip to conform to the shape of your nose



4. Cover the mask with your hands and exhale. If any air leaks from the edge of the mask, adjust the mask again until no air leaks

All information, recommendations, and advice on the part of R&G Faserverbundwerkstoffe GmbH are published to the best of our knowledge and belief. They are noncommittal and contain neither explicit nor tacit assurance or warranty of particular properties. The values specified for properties are typical figures. Recommendations or advice serve to describe our products and possible applications in a general or exemplary, but not specifically individual manner. In the course of the constant technical advancement and improvement of our products there may be changes to the characteristic values, copy, and diagrams; no specific reference is made to any such change. Owing to our products' wide and highly diverse range of potential applications far beyond any of our attempts to analyse, the customer alone is responsible for examining our products' suitability for the respective processes and purposes and their respective processibility. All and any protective rights and the applicable laws, terms, and conditions must be observed by the buyer or user of our products at their own responsibility. Publication is not a licence and does not intend the violation of any patents.



FITTING INSTRUCTIONS

1. With nosepiece facing away from you, hold the bottom strap in each hand with the nosepiece up.
2. Position the respirator under your chin, pull the bottom strap over your head and position it around the neck, below the ears.
3. Pull the top strap over your head resting it high at the top of your head.
4. Place your fingertips from both hands at the top of the metal nosepiece. Using two hands, mold the nose area to the shape of your nose by pushing inward while moving your fingertips down both sides of the nosepiece.
5. Perform a User Seal Check prior to each wearing. To check the respirator-to-face seal, place both hands completely over the respirator and exhale. Be careful not to disturb the position of the respirator. If air leaks around nose, readjust the nosepiece as described in step 4. If air leaks at the respirator edges, work the straps back along the sides of your head.
6. If you CANNOT achieve a proper seal, DO NOT enter the contaminated area. IMPORTANT* Pinching the nosepiece using one hand may result in improper fit and less effective respirator performance. Use two hands.

WARNINGS



Do not reuse



Do not use if package is damaged



Consult instructions for use



Caution consult accompanying documents



Keep away from direct sunlight



Keep dry

CLEANING INSTRUCTIONS

1. Do not alter, modify, clean or repair this respirator.
2. Dispose of the mask in accordance with local regulations.

STORAGE AND TRANSPORT

1. Store according to the manufacturer's instructions.
2. The storage temperature should be between -20 °C and +40 °C.
3. The storage humidity should be below 60 %.
4. When storing or transporting this respirator, use the original packaging supplied with it.
5. Do not store in direct sunlight.

CLASSIFICATION

Three classes are available: FFP1 / FFP2 / FFP3 (see labelling on the packaging and mask). The inscription "NR" on the folded disposable respirator mask stands for "Not reusable" and must be disposed of after single use.

APPROVALS

These respirators are tested by the SGS institution and verified by Universal Certification and Staunchly Management & System Services Ltd.

These respirators are CE marked in accordance with the requirements of European Regulation (EU) 2016/425. The applicable legislation can be obtained on request by checking the certificate and declaration of conformity.

DESIGNATED BODY

Universal Certification and Surveillance Service Trade co., Ltd.
UNIVERSAL SERTİFİKASYON VE GÖZETİM HİZMETLERİ TİCARET LTD. STİ.
Necip Fazıl Bulvarı Keyap Sitesi E2 Blok No:44/84 Yukarı Dudullu Ümraniye /
İSTANBUL / TÜRKİE

All information, recommendations, and advice on the part of R&G Faserverbundwerkstoffe GmbH are published to the best of our knowledge and belief. They are noncommittal and contain neither explicit nor tacit assurance or warranty of particular properties. The values specified for properties are typical figures. Recommendations or advice serve to describe our products and possible applications in a general or exemplary, but not specifically individual manner. In the course of the constant technical advancement and improvement of our products there may be changes to the characteristic values, copy, and diagrams; no specific reference is made to any such change. Owing to our products' wide and highly diverse range of potential applications far beyond any of our attempts to analyse, the customer alone is responsible for examining our products' suitability for the respective processes and purposes and their respective processibility. All and any protective rights and the applicable laws, terms, and conditions must be observed by the buyer or user of our products at their own responsibility. Publication is not a licence and does not intend the violation of any patents.